
	Do not reuse; single use only. Dispositivo de un solo uso; no reutilizar Não reutilize; uso único apenas.	Monouso Ne pas réutiliser, dispositif à usage unique Zum einmaligen Gebrauch bestimmt. Nicht wieder verwenden
	Do not use if package or product damaged. No utilizar si el envase o el producto está dañado Não use se o pacote ou o produto estiverem danificados.	Non usare se la confezione o il prodotto sono danneggiati Ne pas utiliser si l'emballage ou le produit est endommagé Nicht benutzen, wenn Instrument oder Verpackung beschädigt ist
	Specifications subject to change without notice. Las especificaciones están sujetas a modificación sin previo aviso. Especificações sujeitas à mudança sem notificação prévia.	Specifiche soggette a modifica senza preavviso. Les spécifications sont susceptibles d'être modifiées sans préavis. Änderungen vorbehalten.

English – Instructions for Use

Product Description

The IRIDEX ENT Probe is a sterile, single use, disposable fiber optic delivery device for use in surgical procedures for use in incision, excision, coagulation, and vaporization of soft and fibrous tissue including osseous tissue in the medical specialty of ENT surgery. This product should only be used by physicians trained in the use of surgical lasers who are thoroughly familiar with the operation of the laser system and have knowledge of the applicable laser-tissue interaction.

Indications for Use

The IRIDEX ENT Probe is intended for use in surgical procedures for use in incision, excision, coagulation, and vaporization of soft and fibrous tissue including osseous tissue in the medical specialty of ENT surgery. The probe is cleared for use for the particular indications of the compatible laser system to which it is attached.

For information on using this device with a compatible laser system, refer to the laser system operator manual.

Instructions for Use

1. Make sure laser is "OFF" or in "STANDBY" mode prior to attaching the probe to the laser system.
 2. DO NOT USE probe if pouch has been opened or if there is any visual damage to the pouch.
 3. Remove the probe from the pouch and inspect the fiber optic cable for damage.
 4. Remove the protective cap from the probe's fiber optic connector.
 5. Remove protective cap from the laser system fiber port.
 6. Insert the probe into the laser system fiber port and finger-tighten to attach securely.
 7. The probe will display on the laser system panel if properly attached.
 8. The system is now ready to use.
- Note: The aiming beam should be visible when the system is in "TREAT" mode. If not visible refer to the Troubleshooting section of the laser system operator manual.
9. Refer to laser system operator manual for recommended power settings. Always start with lower power settings and short duration times. Increase slowly as necessary for the desired effect.
 10. When treatment is complete, disconnect the probe from the laser system and dispose in accordance with your facilities' biohazard handling procedures.
 11. IRIDEX recommends the safety guidelines of US national standard ANSI Z136.3-1995 as follows:
 - Use the lowest possible oxygen concentration to support the patient;
 - Use the venturi ventilation technique when possible;
 - Use intravenous anesthetic agents rather than inhalation techniques;
 - Use non-flammable laser-safe endotracheal tubes; and
 - Protect the endotracheal tube cuff with wet cottonoids.

Precautions

- CAUTION – This device must be used with appropriate laser safety filter.
- United States Federal law restricts this device to sale by or on the order of a physician.
- This product should only be used by physicians who are trained in, and familiar with, the proper use and operation of the laser system and have knowledge of the applicable laser-tissue interaction.
- Reuse of this product voids all claims of sterility, performance, safety and effectiveness.
- Keep the fiber optic tip free of finger prints and debris and protect from damage.

Warnings

- Proper eye protection must be worn by all personnel and patients in the treatment area to avoid unintended exposure to laser radiation and unintended damage to the patient, user, or environment. Refer to the laser system operator manual for additional information regarding ocular protection.
- The device must be handled with care to prevent damage to the optical fiber [kinking, bending sharply (diameter less than 15 cm or 6 inches), clamping the fiber, etc.] that may result in unintended exposure to laser radiation and unintended damage to the patient, user, or environment. Discard devices that show signs of damage.
- Do not open sterile package prematurely. Open sterile package just prior to use to maintain sterility of contents.

- Verify that the probe is properly connected to the console. An improper connection may result in unintended exposure to laser radiation and unintended damage to the patient, user, or environment.

Potential Complications

The known complications of laser treatment include a possible complication of swelling (edema) in the area treated by the laser. As with any invasive surgical procedure, potential risks of infection and fever exist.

Service and Technical Support

Contact your local IRIDEX Technical Support representative or corporate headquarters should you require assistance.

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