


Fovea-Friendly™ MicroPulse® Laser Therapy



MicroPulse = Increased Retinal Sensitivity without Damage



Prospective, Masked, Randomized Clinical Trial¹

- 62 eyes (50 patients)
- Untreated, center-involving CSME
- Randomized to mETDRS or 810 nm MicroPulse


1 Year Results

- MicroPulse was as effective as mETDRS in
 - stabilizing VA
 - reducing macular edema
- With added benefits of
 - no tissue damage detectable at any time point postoperatively
 - significant improvement in retinal sensitivity

1. Vujosevic S, Bottega E, Casciano M, Pilotto E, Convento E, Midena E. *Retina* 2010



MicroPulse Stimulates Biological Factors



MicroPulse laser treatment, produces a stress response and induces beneficial intracellular biological factors that are primarily anti-angiogenic and restorative without tissue damage as seen in CW.

PEDF - plays a role in inhibiting neovascularization by its anti-angiogenic activity

TSP1 - one of the most potent anti-angiogenic factors

SDF1 - plays a key role in recruitment of bone marrow-derived reparative cells

β -actin - protein that is involved in cell motility, structure and integrity

1. Ogata N, Tombran-Tink J, Jo N, Mrazek D, Matsumura M: Upregulation of Pigment Epithelium-Derived Factor after Laser Photocoagulation. *Am J Ophthalmol* 2001;132(3):427-9.
2. Binz N, Graham CE, Simpson K, Lai YK, Shen WY, Lai CM, Speed TP, Rakoczy PE: Long-Term Effect of Therapeutic Laser Photocoagulation on Gene Expression in the Eye. *FASEB J* 2006;20(2):383-5.
3. Yu AK, Merrill KD, Truong SN, Forward KM, Morse LS, Telander DG: The Comparative Histologic Effects of Subthreshold 532- and 810-Nm Diode Micropulse Laser on the Retina. *Investigative Ophthalmology & Visual Science* 2013;54(3):2216-2224.



Considerations for Incorporating MicroPulse

- Safety
 - Fovea-friendly, no tissue damage, repeatable
 - Absence of laser scars minimizes vision loss over time
- Efficacy
 - Demonstrated through clinical studies and practical experience, durable therapy
- Efficiency
 - Quick and easy treatment
 - Faster than conventional focal or grid laser
- Economics
 - Patient **\$↓**
 - Healthcare System **\$↓**
 - Expense to Practice **\$↓**
 - Dual-laser Platform **\$↓**
 - Practice Revenue **\$↑**



MicroPulse Efficacy & Safety: 810 nm

- Long-term retrospective review: 274 consecutive eyes with macular edema due to DME or BRVO were treated with MicroPulse high density laser treatment using various duty cycles (DC) and followed for up to 10 years. 252 eyes met inclusion criteria.


- Results:

Frequency of laser-induced retinal damage:


- Eyes treated with 10-15% DC 8% (7/84)
- Eyes treated with 5% DC 0% (0/168)

- **5% DC treated eyes showed no detectable retinal damage using infrared, red-free or FAF photos; FA, ICGA; or SD-OCT at 12 months**

MicroPulse Treatment
for CSDME



Pre-treatment



3 Months Post-treatment

Images compliments of
Dr. Jeffrey Luttrull



MicroPulse Efficacy & Safety: 810 nm & 577 nm

- Both 810 nm and 577 nm subvisible MicroPulse laser with 5% duty cycle and fixed power parameters appear to be safe in center involving DME.
- At 6 months, 60 eyes (43 patients) treated with 810 nm and 577 nm showed:
 - No difference in macular volume
 - No signs of inner or outer retinal and choroidal damage
 - No changes shown on FAF or MP1
 - No absolute scotoma
 - Fixation was always central and stable in all patients

Treatment Parameters		
Wavelength	810 nm	577 nm
Eyes	31	29
Spot Size	125 μ m	100 μ m
Power	750 mW	250 mW
Duration	200 ms	200 ms
Duty Cycle	5%	5%

"[MicroPulse] should be introduced into the routine treatment protocols of DME as primary treatment in mild DME (<400 μ m) or combined with intravitreal injections (both anti VEGF or corticosteroids) in moderate to severe DME."




Vujosevic S, Martini F, Convento E, Longhin E, Pilotto E, Midena E: Morphologic and Functional Effects of Diode (810nm) and Yellow (577nm) Subthreshold Micropulse Laser in Center-Involving Diabetic Macular Edema. *ARVO Meeting Abstracts* 2013;54(6):2380.



What is MicroPulse Technology?

MicroPulse technology finely controls thermal elevation by “chopping” a continuous-wave (CW) beam into an envelope of repetitive short pulses.




How Does MicroPulse Work?


Repetitive short pulses permit tissue to cool between pulses and reduce thermal buildup.




CW Laser Exposure
100% Duty Cycle (DC)



MicroPulse High DC (15%)



MicroPulse Medium DC (10%)




MicroPulse Low DC (5%)



MicroPulse Low Intensity/High Density Application

Low-intensity MicroPulse exposures avoid thermal retinal injury. Therefore, **high-density** (confluent) coverage of the diseased retina is needed to maximize clinical effectiveness of MicroPulse Laser Therapy.¹⁻⁸

Continuous-wave Laser




High-intensity argon



Low-intensity argon



Pattern Scanning




Area of retina damaged by laser

Area of retina affected by laser but not destroyed; able to contribute to the therapeutic effects of laser treatment


MicroPulse Laser



Low-intensity/
low-density




Low-intensity/
high density



Area of retina affected by laser but not destroyed; able to contribute to the therapeutic effects of laser treatment




Evolution of Subthreshold Photocoagulation




Barely Visible CW Scanning \neq MicroPulse

Continuous-Wave Barely Visible Scanning



Fluid-like spikes adjacent to ruptured RPE (arrows)
Immediately after CW and pattern scanning laser, a
hyper-reflective band appeared at the laser sites.

IRIDEX MicroPulse



100% of tissue spared Retinal morphology did not change at
any time during the observation period after MicroPulse.


Inagaki K, et al. Spectral-Domain Optical Coherence Tomography Imaging of Retinal Changes after Conventional Multicolor Laser, Subthreshold Micropulse Diode Laser, or Pattern Scanning Laser Therapy in Japanese with Macular Edema. *Retina* 2012;32(8):1592-160



Subvisible CW Scanning \neq MicroPulse

Subthreshold CW Laser

561 nm; 200 μ m spot; 50 mW; 10 ms



Subthreshold CW Laser

532 nm; 200 μ m spot; 50 mW; 30 ms



Photos compliments of Sam Mansour, MD



MicroPulse Clinical Results - Lavinsky

Comparison of mETDRS vs. Low Density vs. High Density Protocols for DME

- A prospective, double-masked, controlled clinical trial on 123 eyes with DME
- Compared three dosing protocols and followed patients for a minimum of 1 yr
- Results:

	Modified ETDRS	MicroPulse High Density	MicroPulse Low Density
Treatment Intensity	Mild	Low	Low
Treatment Density	Low	High	Low
OCT-CMT (Δ)	-126 μm	-154 μm	-32 μm
BCVA (Δ letters)	+4	+12*	-1
Gain ≥ 15 letters	23%	48%*	5%

*Indicates significant improvement compared to mETDRS ($P < 0.05$)



Lavinsky D, Cardillo JA, Melo LA, Jr., Dare A, Farah ME, Belfort R Jr. Invest Ophthalmol Vis Sci 2011





Clinical Success Using Confluent Spacing

Confluent, High-Density Laser Patterns For MicroPulse® Protocols

MicroPulse laser therapy has shown clinical success using confluent spacing.¹⁻²

TxCell™ Scanning Laser Delivery System offers confluent, high-density applications in a wide selection of patterns.





1. Luttrull JK, Sramek C, Palanker D, Spink CJ, Musch DC. *Retina* 2012;32(2):375-86

2. Lavinsky D, Cardillo JA, Melo LA, Jr., Dare A, Farah ME, Belfort R Jr. *Invest Ophthalmol Vis Sci* 2011; 52 (7): 4314-23





TxCell-Guided MicroPulse for Center-Involving ME due to BRVO

Patient: 79-year-old female with BRVO and AMD.



Pre-MicroPulse: VA 20/30
Central retinal thickness 263 μm
Extrafoveal thickness 321 μm
Juxtafoveal thickness 352 μm




15 Weeks Post-MicroPulse: VA 20/25
Central retinal thickness 259 μm
Extrafoveal thickness 312 μm
Juxtafoveal thickness 329 μm

[Read Case Report by Dr. David Gossage](#)




TxCell-Guided MicroPulse for Central Serous Retinopathy


Patient: 54-year-old male with history of bullous CSR and a large persistent neurosensory detachment for 10 months after initial presentation.



Pre-MicroPulse: VA 20/200
CMT 640 μm



4 mos Post-MicroPulse: VA 20/25
CMT 204 μm with complete
resolution of subretinal fluid




TxCell-guided MicroPulse: three 7x7
treatment grids placed confluent
over all areas of leakage or fluid,
including in the fovea.

[Read Case Report by Dr. Gennady Landa](#)



TxCell-Guided MicroPulse for DME

Patient: 79-year-old female with insulin-dependent diabetes mellitus and a history of nonproliferative diabetic retinopathy OU.



Pre-MicroPulse: VA 20/60
CST 340 μm




1 mo Post-MicroPulse: VA 20/40*
CST 280 μm *Baseline for patient

[Read Case Report by Dr. Johnny Tang](#)




Center-Involving DME Refractive to Avastin and Ozurdex

Patient: 63-year-old male with very aggressive refractive DME that would not respond to anti-VEGF or steroid



Pre-MicroPulse: VA 20/200, CRT 438 μm




3 Months Post-MicroPulse: VA 20/60, CRT 270 μm

[Read Case Report by Dr. Aaron Appiah](#)




Macular Edema Associated with CRVO


Patient: 64-year-old male with history of systemic hypertension



Pretreatment: VA 20/150, CRT 870 μm .
Clinical exam revealed prominent cystoid macular edema.



6 weeks post third anti-VEGF treatment.
Pre-MicroPulse: VA 20/70-2, CRT 584 μm .
Recurrent macular edema noted.




Approximately 5 months Post-MicroPulse:
VA 20/40+2, CRT 261 μm . No macular edema observed on clinical exam.

[Read Case Report by Dr. Patrick Caskey](#)



How would you treat these patients with good vision?

All patients have 20/20 to 20/40 VA.



Intrafoveal cysts without retinal thickening.

Intrafoveal thickening with minimal central foveal thickening

Isolated central foveal cyst

Diffuse macular thickening including the fovea


Transfoveal MicroPulse for DME for Patients with Good Vision

Luttrull. RETINA

SAFETY OF TRANSFOVEAL SUBTHRESHOLD DIODE MICROPULSE LASER FOR FOVEA-INVOLVING DIABETIC MACULAR EDEMA IN EYES WITH GOOD VISUAL ACUITY

JEFFREY K. LUTTRULL, MD,* STEPHEN H. SINCLAIR, MD† **RETINA** 34:2010–2020, 2014

Before
4 – 7 mos
Post
MicroPulse



Intrafoveal cysts without
retinal thickening

Intrafoveal thickening with
minimal central foveal thickening

Isolated central foveal cyst


Diffuse macular thickening
including the fovea



Laser Still Plays a Critical Role in the Treatment of DME

In Protocol 1, patients with deferred* laser did best

(*deferred= waiting 6 months before treating with laser)



Considerations for Incorporating MicroPulse

• Safety

- Fovea-friendly, no tissue damage, repeatable
- Absence of laser scars minimizes vision loss over time

• Efficacy

- Demonstrated through clinical studies and practical experience, durable therapy

• Efficiency

- Quick and easy treatment
- Faster than conventional focal or grid laser

• Economics


- Patient **\$↓**
- Healthcare System **\$↓**
- Expense to Practice **\$↓**
- Dual-laser Platform **\$↓**
- Practice Revenue **\$↑**



Repeatable MicroPulse[®] Laser Trabeculoplasty




MicroPulse® Laser Trabeculoplasty (MLT) for the Treatment of Glaucoma



Trabecular meshwork after ALT

Continuous-wave laser exposures can cause high thermal rise resulting in tissue damage




Trabecular meshwork after MLT

Meshwork remains intact without the signs of tissue damage while still as effective as ALT & SLT.¹

1. Fudemberg SJ, Myers JS, Katz LJ, et al: Trabecular Meshwork Tissue Examination with Scanning Electron Microscopy: A Comparison of MicroPulse Diode Laser (MLT), Selective Laser (SLT), and Argon Laser (ALT) Trabeculoplasty in Human Cadaver Tissue. *Invest. Ophthalmol. Vis. Sci.* 2008;49(5):1236-.



Evidence of Dose Response




Multi-Center 1000 mW MLT

Pre-op ≥ 21 mmHg

~ 6 mmHg drop
~ 25% IOP reduction

1 Year IOP Results



Meds

3.0 → 2.8 (NS)

# of pts	Baseline	1 month	3-4 months	6 months	1 year
# of pts	40	34	34	13	10

Ike Ahmed, MD, FRCSC, University of Toronto, Toronto, Canada
David Gossage, DO, FAOC, FAAO, Gossage Eye Institute, Hillsdale, MI
Nathan Radcliffe, MD, Weill Cornell Medical College, New York, NY
Steven Vold, MD, Vold Vision, LLC Fayetteville, AR




Multi-Center 1000 mW MLT

Pre-op ≥ 18 mmHg



1 Year IOP Results



~ 4.5 mmHg drop
~ 20% IOP reduction

Meds

2.9 → 2.3 (NS)

# of pts	Baseline	1 month	3-4 months	6 months	1 year
# of pts	66	54	52	23	14



Ike Ahmed, MD, FRCSC, University of Toronto, Toronto, Canada
David Gossage, DO, FAOC, FAAO, Gossage Eye Institute, Hillsdale, MI
Nathan Radcliffe, MD, Weill Cornell Medical College, New York, NY
Steven Vold, MD, Vold Vision, LLC Fayetteville, AR



MLT / SLT Comparison

	MicroPulse® Laser Trabeculoplasty (MLT)	Selective Laser Trabeculoplasty (SLT)
Wavelength	532 nm and 577 nm	532 nm
Mechanism	Thermally effects - not destroys - pigmented trabecular meshwork cells without thermal or collateral damage	Selective destruction of pigmented trabecular meshwork cells without thermal or collateral damage
Learning Curve	Easy	Easy
Repeatable	Yes	Yes
Visible signs of treatment intra-or post-operative	No	Yes
Inflammation	No	Yes
Spot Size	300 µm (smaller spot to access narrow angles)	400 µm
Complications	Minimal to none	Post-op IOP spikes are possible
Functionality of laser system	Continuous-wave and MicroPulse treatment for glaucoma and retinal disorders	SLT
Parameter Control	Power, ON/OFF time, number and rep rate of pulses	Pulse energy



MLT / SLT Comparison

Micropulse Laser Trabeculoplasty After Previous Laser Trabeculoplasty

BY TAK YEE TANIA TAI, MD

CASE PRESENTATION

A 67-year-old man was referred to me for advanced primary open-angle glaucoma. The patient's visual acuity was 20/25 in the right eye and count fingers at 3 feet in the left eye. He had very mild cataracts in both eyes. A Humphrey 10-2 visual field test (Carl Zeiss Meditec) showed severe constriction that was greater in the left eye. Advanced cupping of the optic disc was present in both eyes, and the IOP was 20 mm Hg in each eye (Figures 1 and 2).

The patient noted that he had been using timolol-brimonidine (Combigan; Allergan) and travoprost (Travatan; Alcon) in both eyes for an extended period of time. Considering the advanced nerve damage, I felt the IOP needed to be lowered further. The patient lived

in Jamaica and traveled frequently, so I first attempted to maximize his medical regimen as much as possible. I started him on methazolamide 50 mg once daily (he was unable to tolerate more frequent dosing), but the IOP in both eyes remained in the high teens.

In February 2013, I performed selective laser trabeculoplasty (SLT) on the patient's left eye. I treated 270° with 75 spots, ranging from 0.5 to 0.7 mJ per spot. I decided against a 360° SLT treatment due to the potentially higher risk of an IOP spike after this procedure with a greater area of laser application.¹ The patient responded well, and the IOP decreased to 13 mm Hg in the left eye. I treated the right eye with the same protocol in April 2013, after which the IOP in both eyes measured between 12 and 13 mm Hg. Because of the patient's advanced optic nerve damage

Case Report

- Patient IOP rebound post SLT
- Physician concern with SLT repeat treatment because of potential IOP spikes
- MLT performed because no history of inflammation post treatment
- Patient's IOP dropped from 19 to 13 mm Hg and stable for 6 months

"For patients such as this one, with advanced disease and IOP spikes and for whom pharmaceutical treatments and previous laser treatments have failed, MLT is a viable option"

[Read Case Report by Dr. Tania Tai](#)



MicroPulse Technology available in Multi-Functional Laser Systems



IQ 577™



IQ 532™



IQ 810™

- Fovea-friendly™ MicroPulse* Laser Therapy¹ for retinal disorders
- Repeatable MicroPulse Laser Trabeculoplasty for glaucoma therapy
- Conventional photocoagulation
- TxCell™ Scanning Laser Delivery Device*: Multi-spot pattern scanning for efficient retinal photocoagulation
- FDA and CE clearance for both conventional and MicroPulse Laser Therapies

*MicroPulse and TxCell are optional with the IQ 577 and IQ 532 lasers
1. Bhagat N, Zarbin M, Mansour S, Chong V, and Cardillo JA. Fovea-friendly MicroPulse Laser. Supplement to Retina Today, May/June 2012



New Standard of Care

