IRIDEX MicroPulse™ P3 Device Instrumentation

English

Instructions for Use

INDICATIONS:

Indications for this device include, but are not limited to transscleral photocoagulation for the treatment of primary open-angle glaucoma, closed-angle glaucoma, and refractory glaucoma using MicroPulse™.

DIRECTIONS FOR USE:

Remove probe from package and unwind with care. This product contains a glass optical fiber that can be damaged with improper handling. Connect probe to a compatible IRIDEX laser console. Consult your operator manual for additional instructions, contraindications, warnings, and cautions.

Anesthesia

Typically an anesthesia block is used during these procedures. Keep the tip and the eye surface moist throughout procedure.

Placement (Figure 1)

Place the shorter radius of the distal tip of the device concentric to the limbus, and hold it approximately perpendicular, see Figure 1. Note, the device will align the internal fiber-optic tip at approximately 3mm posterior to the limbus.

Fig. 1: Placement

Application

- 1) The following method of application is intended to provide guidance only for treatment settings which are not prescriptive for any condition. The operative needs of each patient should be individually evaluated based on the specific indication, treatment location, and patient-specific characteristics. If uncertain of expected clinical response, always start with conservative settings and increase laser power and/or duration settings in small steps.
- 2) To treat the superior quadrant, apply firm pressure and move the device in a continuous sliding arc motion from the 9:30 to 2:30 clock positions. Apply laser continuously, "PAINT", for 4-7 complete passes between the clock positions, for a total of 50 seconds. Based on recommendations of experienced clinicians: use MicroPulse mode at 31.3% (0.5 ms on 1.1 ms off).¹
- 3) If device snags on conjunctiva, momentarily stop laser treatment to release caught conjunctiva and reapply the device at same spot to continue treatment.
- 4) If a pre-existing trabeculectomy, tube, shunt, or area of scleral thinning exists, reduce the clock hours to avoid that area. Still deliver a total of 50s over the reduced clock hour range.
- 5) To treat the inferior quadrant, repeat process for a total of 50s over the 3:30 to 8:30 clock hours.

Device Fiber Cleanliness

Keep the device tip clean to minimize the risk of burns to the ocular surface. If the tip accumulates debris or "charring" during the procedure, clean it gently with an alcohol swab. If the "charring" or discoloration on the tip cannot be removed by gentle cleaning, discard the device. Scleral burns are not typical and may indicate contamination at the device tip. If a scleral burn occurs, discontinue use and replace the device immediately. The device is a Single-Use Product.

CAUTION:

Keep the device tip and the eye surface moist throughout treatment. Do not treat over sites of previous trabeculectomies, scleral thinning, nor directly over tubes and shunts.

WARNING:

Excessive treatment power may result in ocular surface burns or ciliary body hemorrhage. Contamination of the fiber optic tip by blood or tissue char may result in ocular surface burns. Excessive energy may cause equatorial burns. Heavy perilimbal conjunctival pigmentation may result in local absorption and burns; therefore, avoid areas of heavy perilimbal pigmentation. This device must be used with appropriate laser safety filter or eyewear. Never look directly into the laser light source or at laser light scattered from reflective surfaces.

Do not open sterile package prematurely. Open sterile package just prior to use to maintain sterility of contents.

Inspect packaging prior to use: DO NOT USE IF THE PACKAGE IS DAMAGED OR IF THE STERILE BARRIER IS COMPROMISED.

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IRIDEX Corporation