



Operating Microscope Adapter (OMA)

Features and Benefits

- Offers precise targeting and the therapeutic capability of retinal photocoagulation
- Spot sizes include 0.3, 0.5, 0.8, 1.2, and 2.0 mm
- Rugged construction makes it ideal for an operating room environment
- Compatible with a variety of operating microscopes

Microscope Mfg.	Microscope Model	OMA (IQ 810™)	OMA (SL/SLx)	Comments
Zeiss	Lumera	60429	13411	All Zeiss models require Zeiss dovetail eye safety filter #10378
	i			
	700			
	T			
	S7			
	Visu	N/A	N/A	
	140 [^]			
	150 [^]			
	160 [^]			
	200 [^]			
	210 [^]	60429	13411	One assistant beam path will be obstructed.
	CS w/variospot	60436	11674	
	CS w/retrolux-3	60434	11890	
	MDO, MDU	60433	11673	
OPMI Series 1, 19	60431	11688		
OPMI Series 6	60435	11672		
OPMI 1 FR pro	60429	13411	Requires Zeiss accessory #303360-9910	
Wild/Leica	M840, M841, M844	60699	14308	Requires Wild/Leica dovetail eye safety filter #14080 or Qty 2 of #11175. One assistant beam path will be obstructed.
	M500, M501, M650, M651, M690, M691	60430	11671	Requires Wild/Leica dovetail eye safety filter #11175
Weck	OM-1205	60432	11818	Requires Zeiss dovetail eye safety filter #10378

[^]Discontinuation confirmed by manufacturer. N/A = Product not available.

Note: All operating microscopes must be equipped with a 175 mm objective lens which can be purchased from the microscope dealer (not from IRIDEX).

Specifications are subject to change without notice. IRIDEX and the IRIDEX logo are registered trademarks.

Products are covered by one or more of the following U.S. patents: 5,511,085; 5,982,789; 6,327,291; 6,540,391; 6,733,490; 7,766,904; 7,771,417; 7,909,816; 5,997,498; 6,073,759; 6,092,898; 6,217,594; 6,494,314; 6,585,679; 6,726,666; 6,800,076; 6,866,142; 7,537,593; 8,177,777; 7,837,83; 6,953,049; 6,766,904; 7,771,417; 7,909,816; 5,997,498; 6,073,759; 6,092,898; 6,217,594; CA 2286002; JP 449444; JP 4570696; JP 4819754; JP 5123973; JP 5133069. Other U.S. and international patents pending.

Clinical references available upon request.

EC REP Emergo Europe
Molenstraat 15, 2513 BH, The Hague, The Netherlands, Tel.: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299

